


# PRODUCT DATA SHEET

Revision 1	96 Well Cell Culture Plate, White, F-Bottom	
Validfor Item-No.:	1014100/1014104/1014160/1014164/1014102/1014106/ 1014162/1014166	

1.	Description / Specification	
1.1	Description	Cell Culture Plate, 96 well, solid white flat bottom , with vents, alphanumeric well coding, single position lid with condensation rings, sterile 1014100 96 well, Cell culture plate, with white lid, TC treated, sterile 1014104 96 well, Cell culture plate, with white lid, TC treated, sterile 1014160 96 well, Cell culture plate, with clear lid, TC treated, sterile 1014164 96 well, Cell culture plate, with clear lid, TC treated, sterile 1014102 96 well, Cell culture plate, with white lid, TC treated, sterile 1014106 96 well, Cell culture plate, with white lid, TC treated, sterile 1014162 96 well, Cell culture plate, with clear lid, TC treated, sterile 1014166 96 well, Cell culture plate, with clear lid, TC treated, sterile
1.2	Dimensions	Plate : see Customer Drawing 1014100-1014164:Individually Wrapped in Peelable 1014102-1014166:Individually paper-plastic wrapped
1.3	Volume	Max volume: 0.36 ml (mathematical calculated) Working volume: 0.07-0.2 ml Growth area / well: 0.32 cm <sup>2</sup>
1.4	Material / Resin	Plate : PS (Polystyrene) Lid: PS (Polystyrene) The materials for manufacturing are <i>USP Class VI</i>
1.5	Colour	Plate :black Lid: clear/ black
1.6	Sterilization	SAL 10 <sup>-6</sup>
1.7	Quality Control	- Raw Material-Control: physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Expiry date and Lot-No. printed on the outer box

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	1
3.2	Pieces / Box	50/100
3.3	Lot-No.	DDD YY XXX (day, year, product No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
4.1	Laboratory use only

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 07.05.2024	Date 14.05.2024	Date 31.05.2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.